
































































Straumann® Dental Implant System System Overview <b>TISSUE LEVEL SURGERY</b> MRK.034.Rev.04		Ø 3.3 mm				Ø 4.1 mm				Ø 4.8 mm											
Titanium SLA®		Standard Plus (SP) RN Ø 4.8 mm  Ø 3.3 mm				Standard Plus (SP) RN Ø 4.8 mm  Ø 4.1 mm				Standard Plus (SP) RN Ø 4.8 mm  Ø 4.8 mm		Standard Plus (SP) WN Ø 6.5 mm  Ø 4.8 mm									
Roxolid® SLActive®		4 mm – 6 mm – 8 mm 043.2605 10 mm 043.2615 12 mm 043.2625 14 mm 043.2635		Roxolid® SLA® only		4 mm – 6 mm 043.9255 8 mm 043.9265 10 mm 043.9275 12 mm 043.9285 14 mm 043.9295				4 mm – 6 mm 043.9395 8 mm 043.9405 10 mm 043.9415 12 mm 043.9425 14 mm 043.9435		4 mm – 6 mm 043.9485 8 mm 043.9495 10 mm 043.9505 12 mm 043.9515 14 mm –									
Round Burs		 Round Bur Ø 2.3 mm 044.003 →  Round Bur Ø 3.1 mm 044.004		 Round Bur Ø 2.3 mm 044.003 →  Round Bur Ø 3.1 mm 044.004		 Round Bur Ø 2.3 mm 044.003 →  Round Bur Ø 3.1 mm 044.004				 Round Bur Ø 2.3 mm 044.003 →  Round Bur Ø 3.1 mm 044.004											
Pilot and Twist Drills		 Pilot Drill Ø 2.2 mm 044.783 short 044.785 long →  Twist Drill Pro Ø 2.8 mm 044.787 short 044.789 long		 Pilot Drill Ø 2.2 mm 044.783 short 044.785 long →  Twist Drill Ø 2.8 mm 044.787 short 044.789 long →  Twist Drill PRO Ø 3.5 mm 044.791 short 044.793 long		 Pilot Drill Ø 2.2 mm 044.783 short 044.785 long →  Twist Drill PRO Ø 2.8 mm 044.787 short 044.789 long →  Twist Drill PRO Ø 3.5 mm 044.791 short 044.793 long →  Twist Drill PRO Ø 4.2 mm 044.795 short 044.797 long				 Pilot Drill Ø 2.2 mm 044.783 short 044.785 long →  Twist Drill PRO Ø 2.8 mm 044.787 short 044.789 long →  Twist Drill PRO Ø 3.5 mm 044.791 short 044.793 long →  Twist Drill PRO Ø 4.2 mm 044.795 short 044.797 long											
Alignment Pins Depth Gauges		 Alignment Pin Ø 2.2 mm 046.799 →  Depth Gauge Ø 22.8 mm 046.800		 Alignment Pin Ø 2.2 mm 046.799 →  Depth Gauge Ø 2.8 mm 046.800 →  Depth Gauge Ø 3.5 mm 046.802		 Alignment Pin Ø 2.2 mm 046.799 →  Depth Gauge Ø 2.8 mm 046.800 →  Depth Gauge Ø 3.5 mm 046.802 →  Depth Gauge Ø 4.2 mm 046.804				 Alignment Pin Ø 2.2 mm 046.799 →  Depth Gauge Ø 22.8 mm 046.800 →  Depth Gauge Ø 3.5 mm 046.802 →  Depth Gauge Ø 4.2 mm 046.804											
Profile Drills		 SP Profile Drill 034.332		 SP Profile Drill 034.333		 SP Profile Drill 044.084				 SP Tap Ø 4.8 mm 034.353		 BL Tap Ø 4.1 mm 034.350 רק לשתול באורך 4 ס"מ									
Taps		 SP Tap Ø 3.3 mm 034.351		 SP Tap Ø 4.1 mm 034.352		 BL Tap Ø 4.1 mm 034.349 רק לשתול באורך 4 ס"מ				 SP Tap Ø 4.8 mm 034.353		 BL Tap Ø 4.1 mm 034.350 רק לשתול באורך 4 ס"מ									
Closure Caps Healing Caps		RN				RN				RN				WN							
		Closure Cap - Flat  048.3715				Healing Abutments  Ø 6 mm H 1.5 mm 048.3735  Ø 5.5 mm H 2 mm 048.0335  Ø 5.5 mm H 3 mm 048.0345  Ø 5.5 mm H 4.5 mm 048.0375								Closure Caps  H 0 mm 048.3755  H 2 mm 048.0385  H 3 mm 048.0395  H 4.5 mm 048.0535				Healing Abutments			



**Impression Posts**  
**Impression Caps**  
**Implant Analogs**



**Scanbodies and Repositionable Analogs**



**Temporary Restorations**







**Prosthetic Restoration** שיקום חלקי

**Screw-retained** מבנה מוברג


**Cement-retained** מבנה מודבק

**Screw-retained** מבנה מוברג

**Cement-retained** מבנה מודבק







**Titanium Abutments**



synOcta® Abutment			
AH	◄	REF	
5.5mm	0°	048.605	
6.7mm	15°A	048.610	
6.7mm	15°B	048.611	
6.7mm	20°A	048.615	
6.7mm	20°B	048.616	



**Plastic Copings**



**CoCr Copings**  
**CoCr Abutments**



**Auxiliary Parts**  
**Screws**



**Kit TL Modular** קיט מודולרי

Adapter for Ratchet 046.460, 046.461  
Adapter for Handpiece 046.470, 046.471  
Drill Extender 040.563  
BL/NNC Tap 034.348  
BL/NNC Profile Drill 034.327  
S/SP Tap 034.351  
SP Profile Drill 034.332  
S/SP Tap 034.352  
SP Profile Drill 034.333  
S/SP Tap 034.353  
SP Profile Drill 044.084  
BL Tap 034.349  
BL Tap 034.350

Ø 3.1, Ø 2.3, Ø 1.6, Ø 2.2, Ø 2.8, Ø 3.5, Ø 4.2, long, short

Round Bur 044.004, 044.003  
Needle Drill 026.0056  
Pilot Drill 044.783, 044.785  
Alignment Pin 046.799  
Twist Drill PRO 044.787, 044.789  
Depth Gauge 046.800  
Twist Drill PRO 044.791, 044.793  
Depth Gauge 046.802  
Depth Gauge 046.804